



ParaFishControl



New tools for the diagnosis of the amoebic gill disease (AGD)

ParaFishControl Workshop “North European Fish Parasite Management Strategies in Aquaculture Farms”

Berlin, 9th October 2019

Ariadna Sitjà-Bobadilla O.B.



Irene Cano, Brian Mulhearn & Richard Paley



International Centre of Excellence for Aquatic Animal Health, Weymouth Laboratory, Weymouth, Dorset, DT4 8UB, United Kingdom



Table of Contents

1. Challenge and impact in the industry
2. Our approach – Our team
3. Our proposed solution
4. Our solution vs currently employed solutions
5. Expected benefits for the industry
6. Status and next steps
7. Conclusions

Challenge and Impact



- Amoebic gill disease (AGD) is caused by *Neoparamoeba perurans*
 - AGD **clinical signs**: multifocal white mucoid spots and **patches on the gill surface**
 - A “**gill score**” scale, internationally used:
 - 0 = no lesions
 - 5 = extensive lesions is:
- The **gill score** is **presumptive**




Photo from Journal of Fish Diseases 2011, 34, 411–432

- AGD causes economic **losses to salmonid industry** worldwide.
- **Costs** estimated: at A\$20m annual (Australia), £30m (UK/EIRE, 2011-2012)
- **Treatments** (e.g. freshwater, H₂O₂, peracetic acid..) are **expensive!!!**
- “**Handling**” of the disease: **by gill scoring and early treatment**
- **An on-site confirmatory test** of the presence of *N. Perurans* is desirable before treatments are conducted **for a cost-effective treatment**

Our approach and our team

The proposed solution:

- To develop a **point of care (POC) test** for the **on-site confirmation of *N. perurans*** on non lethal gill-swabs.
- To compare **POC DNA extraction methods** for field application.

Team involved in the development:



Test design


Richard Paley
Irene Cano

Experimental challenges and sample collection

Ava Wayne
Claire Joiner

Sample analysis

Brain Mulhearn



POC DNA extraction methods

Robin McCullough
Irene Cano


Amoebae *in vitro* culture

All of them!


Our proposed solution (I)




An isothermal amplification (**LAMP**) **assay** using portable **weather proof** real-time fluorescence equipment for the **on-site detection** of the 18S rRNA gene of *N. perurans* and the host EF1a gene (internal control) in **non-lethal** Atlantic salmon **gill swabs**.




AGD LAMP assay optimization, specificity and analytical sensitivity



Optimal isothermal amplification at **62.9 °C**



AGD LAMP assay **highly specific!**



Detects **10³ copies** parasite 18S rRNA gene-35 min




Our proposed solution (II)

DNA extraction methods for field application

Five “fast and dirty” DNA extraction methods were evaluated:

1. Sodium hydroxide (NaOH)
2. QuickExtract™ DNA extraction solution
3. Potassium hydroxide (KOH)
4. KAPA-express-extract
5. Buccalyse DNA release kit


and compared with the reference laboratory method (EZ1 Robot)



- ✓ **QuickExtract** faster POC DNA extraction method
- ✓ **Compatible** with **Taqman qPCR** chemistry - allows for **laboratory validation** of field LAMP results

Our proposed solution (III)

N. perurans detection in non-lethal gill swabs by LAMP




Bars: Atlantic salmon gill swabs analysed
Line: percentage of LAMP-positive tests per gill score

Expected benefits for the industry ParaFishControl

- **On-site discrimination** of the *N. perurans* from **non-lethal gill swabs** in clinically infected Atlantic salmon
- In parallel with gill scores **can help to cost-effectively plan treatments**
- **Easy** assay to perform (low technical requirements)
- **Fast: Results under 40 min** (5 min for DNA extraction and 35 min to detect 10^3 copies of the targeted gene)
- Gross estimated **cost of £2.5 per sample** (including the cost for DNA extraction and the isothermal amplification)

Current status and next steps

- Test ready for **inter laboratory validation**
- To develop a **commercial kit is possible**
- Possible to develop LAMP assays for **other targets** as new sequences of *N. perurans* genes are available, which might help to improve sensitivity and further reduce test time
- Further **refinement** is required for improved detection of **lower gill scores**



Conclusions

- The developed AGD-LAMP assay performs well (**100% of positives**) for swabs from fish **showing gill scores >3**.
- **Lower gill scores** do **require a bigger sample** size to confirm the presence of *N. perurans* in non-lethal gill swabs.
- Overall, the simplicity, performance and low cost of the present AGD-LAMP assay make this test a **good candidate for the on-site confirmation of *N. perurans* in non-lethal gill swabs** taken from AGD clinically infected Atlantic salmon.
- Further **refinement** is required for improved detection of **lower gill scores**.

Thank You



ParaFishControl

Richard Paley & Irene Cano

Richard.paley@cefas.co.uk

Irene.cano@cefas.co.uk

 @parafishcontrol.eu

www.parafishcontrol.eu



This project has received funding from the European Union's Horizon 2020 research and innovation programme under grant agreement No 634429. This output reflects the views only of the author(s), and the European Union cannot be held responsible for any use which may be made of the information contained therein.